



Email : sales@laboratoryinstrumentindia.com

Product Name :
Current Sensor

Product Code :
LBNY-0004-780059



Description :

Current Sensor

Technical Specification :

Current Sensor :

The Current Sensor determines the current through it by measuring the voltage across the internal 1.00 Ω resistor. Up to 1.5 A can be measured.

Applications:

Study resistance, voltage and capacitance in circuits

Specifications:

Current up to 1.5A

Resolution: 5 mA at 1X gain, 0.5 mA DC



Laboratory Instrument India

