



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Journal Bearing Demonstration

**Product Code :**  
LBNY-0005-14000046



**Description :**

To meet propelling demands of our valuable clients, we are presenting ourselves as dependable manufacturers, exporters and suppliers of a wide assortment Journal Bearing Apparatus. Owing to their quality, offered range is highly appreciated around the globe as it is manufactured as per international standards. Assisted by technically advanced facilities, this is designed by our skilled professionals who execute their work in an efficient and precise manner.

**Features:**

Precisely designed  
Low maintenance  
Highly durable.

**Technical Specification :**

Motor Variable speed Motor  
Control Panel For speed control of motor  
Manometer 16 Tubes  
Bearing Diameter 55 mm  
Journal Diameter 50 mm  
Weights 4 adjustable weights.

---

# Laboratory Instrument India

