



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Field and Detector Coils

**Product Code :**  
LBNY-0004-780060



**Description :**

Field and Detector Coils

**Technical Specification :**

Field and Detector Coils

200 turns of #22 copper wire:

Around 18 cm ID & 22 cm OD

Current up to 2A

500 turns of #22 copper wire

Around 18 cm ID & 22 cm OD

Current up to 2A

Coil has the capability to be connected to the interface, function generator or oscilloscope for input and output data.

Coil detector 400 turns

Coil detector 2000 turns

Set of 2 LED indicators



**Laboratory Instrument India**

