



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Fiber Optic Splicer

**Product Code :**  
LBNY-0002-610009



**Description :**

Fiber Optic Splicer

**Technical Specification :**

Viewing method and display 2 CMOS cameras, 4.3 inch color LCD monitor

Splicing program: 40 groups

Power supply Adaptor, input: AC100-240V(50/60HZ)

Typical splicing time 8 seconds

Terminals USB 2.0 port, for uploading splice results and upgrading software

Loss evaluation: Available

Multi-Mode (MM) G651

Single Mode (SM) G652 & G657

Single Mode (DS/NZDS) Non-Zero Dispersion

Shifted Fiber (G65s) Cladding diameter: 80...150µm

Fiber Optic Splicer, Fiber Optic Splicer Equipments, Fiber Optic Splicer Tools, Fiber Optic Splicer Tool Kits, Fiber Optic Splicer Manufacturers, Fiber Optic Splicer Suppliers from India, China, Kenya

---

# Laboratory Instrument India

