



Email : sales@laboratoryinstrumentindia.com

Product Name :
Heat Pump In Air Conditioning Training Bench

Product Code :
LBNY-0005-14100017



Description :

The Heat Pump in Air Conditioning Training Bench has been designed and manufactured to incorporate all current safety standards.
A Heat Pump in Air Conditioning Training Bench for studying the theoretical and practical aspects of the use of heat pumps in air conditioning systems.
The refrigerant used complies with anti-pollution legislation.

Experiments:

Analysing system behaviour and reactions with the variation of condenser and evaporator air flow
Plotting the cycle in a refrigerant pressure-enthalpy diagram
Studying the operation of a heat pump
Studying the operation of a cycle reversing valve.

Technical Specification :

Obtaining data readings and calculation of:
Heat balance on air circuit at the condenser with evaluation of the heat exchange coefficient
Heat exchange surfaces
Ideal and actual Energy Efficiency Ratio
Ideal and actual coefficient of performance
Heat balance on air circuit at the evaporator
Heat balances for the evaporator, condenser and compressor
Refrigerant mass flow
Volumetric compression efficiency.



Laboratory Instrument India