



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Water Condensing Units Training Plant

**Product Code :**  
LBNY-0005-14100024



**Description :**

A Water Condensing Units Training Plant designed for the study of the different methods used for water condensing in refrigeration units.

The allows students to investigate and study the operation of a pressure operated valve connected to a water network, a dry cooler and a cooling tower.

The Water Condensing Units Training Plant has been designed for the study of the different methods used for water condensing in refrigeration units.

The refrigerant used complies with anti-pollution legislation.

The has been designed and manufactured to incorporate all current safety standards.

**Technical Specification :**

Water Condensing Units Training Plant



**Laboratory Instrument India**