



Email : sales@laboratoryinstrumentindia.com

Product Name :
Quad Temperature

Product Code :
LBNY-0004-780063



Description :

Quad Temperature

Technical Specification :

Quad Temperature

Can be connected to 4 temperature probes

Includes 2 stainless steel probe and at least 3 fast response probes

Measure at least from -30 to 105 degrees celcius

Can be used with ideal gas law apparatus

The Quad Temperature Sensor can connect up to four Temperature Probes and can be used with our Stainless Steel, Fast Response, and Skin/Surface Temperature probes for a wider variety of temperature measurements in the classroom or field.

Features:

Displays °C, K and °F with 0.0025 °C resolution.

Rugged, stainless steel probe with easy-grip plastic handle.

Teflon cover to protect the probe from aggressive chemicals is available separately.

Laboratory Instrument India

