



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Domestic Air Conditioning Assembly Kit

**Product Code :**  
LBNY-0005-14100032



**Description :**

The has been designed and manufactured to incorporate all current safety standards.  
The refrigerant used complies with anti-pollution legislation.  
A Domestic Air Conditioning Assembly Kit for diagnosing problems associated with refrigeration circuits and air conditioning systems.  
The kit comes pre-filled with refrigerant.

**Technical Specification :**

System start-up  
Carrying out refrigerant gas recovery procedures  
Analysing system operating parameters  
Procedures for:  
Filling with refrigerant gas and checking joints  
Draining and cleaning the system.



**Laboratory Instrument India**

