



Email : sales@laboratoryinstrumentindia.com

Product Name :
Domestic Air Conditioning Assembly Kit

Product Code :
LBNY-0005-14100032



Description :

The has been designed and manufactured to incorporate all current safety standards.
The refrigerant used complies with anti-pollution legislation.
A Domestic Air Conditioning Assembly Kit for diagnosing problems associated with refrigeration circuits and air conditioning systems.
The kit comes pre-filled with refrigerant.

Technical Specification :

System start-up
Carrying out refrigerant gas recovery procedures
Analysing system operating parameters
Procedures for:
Filling with refrigerant gas and checking joints
Draining and cleaning the system.



Laboratory Instrument India

