



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hydraulic Servo Universal Testing Machines High Accuracy

**Product Code :**  
LBNY-0005-1420004



**Description :**

**Features:**

- Support for multi-language user-interface extensibility.
- Complies with standard requirements of torsion and compression, as well as a variety of destructive, non-destructive and multi-cycle tests.
- Support for multi-level elastic modulus tests.
- Characteristic testing data can be marked on testing graphs.
- Capable of statistical operations on testing data (average, standard deviation, coefficient of variance).
- Automatic requesting for testing result data, also capable of launching an advanced analysis to inspect the data.
- Dynamic display of load, displacement, deformation and real-time testing curve on the screen.
- Closed loop control of channels such as load, displacement, deformation, etc. Smooth switching in between.
- Rounding off of numerical testing data complies with Rules of Rounding Off for Numerical Values, supports rounding off in sections and custom rounding off.
- As many as 20 layers of sensor linear compensation, elevates the precision of the instruments to the next level.

**Technical Specification :**

Hydraulic Servo Universal Testing Machines High Accuracy

---

# Laboratory Instrument India

