



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Computerized Hydraulic Universal Testing Machine  
Strength Test

**Product Code :**  
LBNY-0005-1420006



**Description :**

Specially designed extended grip jaws are equipped for pre-loading tension on concrete steel strand. This series is also capable of both static tensile, compression tests and quasi-dynamic tests including low-cycle fatigue tests. Apply appropriate grips to perform tension on metal bar, plate and threaded iron. Meanwhile, applicable grips provide capability of tension and compression on metals or concrete.

The mechanical components collaborate with the pioneered modular Test Live computer testing software to ensure the stability, precision and reliability of the whole system. This product features a firm structure as well as steady loads, precise controls and is intuitive to operate. It adopts the leading high-performance Controller, which is usually equipped on high-end material testing machines.

**Technical Specification :**

Computerized Hydraulic Universal Testing Machine Strength Test



**Laboratory Instrument India**