



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Absolute Pressure Sensor

**Product Code :**  
LBNY-0004-780064



**Description :**

Absolute Pressure Sensor

**Technical Specification :**

Absolute Pressure Sensor

This sensor is ideal for the investigation of gas laws. It measures absolute pressure so that in a normal room environment it will measure the atmospheric pressure which is 100KPa. It can therefore measure both a vacuum and higher pressures. When used with the Volume Sensor it is easy to investigate the relationship between pressure and volume.

Pressure range : 0-700 kPa, at least

Provide necessary accessories for air compression, (PVC tube, syringe, etc.)



**Laboratory Instrument India**