



Email : sales@laboratoryinstrumentindia.com

Product Name :
Absolute Pressure Sensor

Product Code :
LBNY-0004-780064



Description :

Absolute Pressure Sensor

Technical Specification :

Absolute Pressure Sensor

This sensor is ideal for the investigation of gas laws. It measures absolute pressure so that in a normal room environment it will measure the atmospheric pressure which is 100kPa. It can therefore measure both a vacuum and higher pressures. When used with the Volume Sensor it is easy to investigate the relationship between pressure and volume.

Pressure range : 0-700 kPa, at least

Provide necessary accessories for air compression, (PVC tube, syringe, etc.)



Laboratory Instrument India