



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Servo Universal Testing Machines With High Accuracy No Interstice

Product Code :
LBNY-0005-1420008



Description :

Features:

Automatic requesting for testing result data, also capable of launching an advanced analysis to inspect the data. As many as 20 layers of sensor linear compensation, elevates the precision of the instruments to the next level. Support saving, instant retrieving and inspection of testing results. Capable of re-analyzing test data and curves, and partial zooming of testing curves at will. Rounding off of numerical testing data complies with Rules of Rounding Off for Numerical Values, supports rounding off in sections and custom rounding off. Characteristic testing data can be marked on testing graphs. Capable of statistical operations on testing data (average, standard deviation, coefficient of variance). Dynamic display of load, displacement, deformation and real-time testing curve on the screen. Closed loop control of channels such as load, displacement, deformation, etc. Smooth switching in between. Support Excel report templates and customized editing of test reports. Integrated with software limit protection, hardware travel limit protection and full scale overload protection. Support for multi-language user-interface extensibility.

Technical Specification :

Hydraulic Servo Universal Testing Machines With High Accuracy No Interstice

Laboratory Instrument India

