




Email : sales@laboratoryinstrumentindia.com

Product Name :
Electro Hydraulic Servo Shear Testing Machine

Product Code :
LBNY-0005-1440003



Description :

Multi-channel coordinative digital load control system makes it more convenient to perform tests. The hydraulic system features a quick withdraw unit to promote efficiency. Electromechanically adjustable test space between upper and lower compression platen. Horizontal shear components are moveable electromechanically and easy to operate. Easy to set testing methods. Administrators may configure in a specific range of settings under related standards, pre-installed with two sets of testing schemas (standard testing schema and custom testing schema). Support for management permissions for various user types (Operator, administrator and engineer). Highly automatic testing process. Capable of general functionalities including automatically aligning with zero, smooth continuous loading and un-loading, automatic sustaining. Capable of monitoring the whole testing process, automatic testing data measuring and processing, automatic testing curve displaying, automatic saving the source data and curves of tests. A wide range of test projects are supported by the testing software (14 categories), perfect support for the latest mechanics performance testing standards on rubber bearing products.

Technical Specification :

Electro Hydraulic Servo Shear Testing Machine

Laboratory Instrument India

