




Email : sales@laboratoryinstrumentindia.com

Product Name :
Microcomputer Control Electro

Product Code :
LBNY-0005-1440004



Description :

Capacity 5000kN
Load measurement range 0.4%100%
Load measurement accuracy 1%
Load resolution 1/500000 of Max capacity, no stepping, no change in full scale
Piston displacement measurement accuracy 0.5%
Piston speed adjustment range 0.140mm/min
Piston travel of vertical cylinder 200mm
Horizontal shear capacity 1000kN
Horizontal piston travel 250mm
Rotation lifting capacity 300kN
Upper/lower compression platen distance 1000mm(or customized)
Compression platen dimensions 820mm820mm
Axial deformation measurement range 010mm
Lateral deformation measurement range 010mm
Horizontal shear deformation measurement range 0150mm
Frame dimensions 490012003000mm3
Hydraulic system dimensions 1040720975mm3
Power supply 11kW/AC380V/50Hz

Technical Specification :

Microcomputer Control Electro



Laboratory Instrument India