



Email : sales@laboratoryinstrumentindia.com

Product Name :
Test Equipment To Measure Weight Length Minus
Tolerance

Product Code :
LBNY-0005-1450003



Description :

This product is well structured, easy to operate, integrated with reliable safety barriers. The measurement data are retrieved and processed by computer software. The Weight Length Minus Tolerance Measuring Instrument is mainly designed to measure the length and weight of steel components.

Technical Specification :

Length measurement range: 2001500mm
Length measurement accuracy: 0.2mm
Weight measurement range: 0.330kg
Weight measurement accuracy: Less than 0.5% of total weight
Sample diameter range: 550mm
Power supply: 0.5kW/AC220V/50Hz
Frame dimensions: 21807201175mm3
Weight: Approx. 750 kg



Laboratory Instrument India

