



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Test Equipment To Measure Weight Length Minus  
Tolerance

**Product Code :**  
LBNY-0005-1450003



**Description :**

This product is well structured, easy to operate, integrated with reliable safety barriers. The measurement data are retrieved and processed by computer software. The Weight Length Minus Tolerance Measuring Instrument is mainly designed to measure the length and weight of steel components.

**Technical Specification :**

Length measurement range: 2001500mm  
Length measurement accuracy: 0.2mm  
Weight measurement range: 0.330kg  
Weight measurement accuracy: Less than 0.5% of total weight  
Sample diameter range: 550mm  
Power supply: 0.5kW/AC220V/50Hz  
Frame dimensions: 21807201175mm3  
Weight: Approx. 750 kg



**Laboratory Instrument India**

