



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Display Electronic Universal Testing Machine

Product Code :
LBNY-0005-1470003



Description :

The test adopt high accuracy and digital speed system, precision decelerate structure and precision drive screw ball to realize the max range of test speed. During the testing there is lower noise and smooth operation. This machine equips alternating servomotor for driving crosshead and imported photoelectric encoder for displacement measure. This machine is accompanied with the digital multi-channel closed loop control system, which apply with the latest measurement and control systems with lots of advanced professional and technical development.

It can achieve the adjustment of crossheads speed when clamp the specimen.

Adopt the USB connector . it can connect the computer to achieve the datas memory, analyze, print and display. Touch button operation, LCD display screen. It includes test methods display screen, test force display screen, test operation and result display screen and curve display screen. Its so convenient and fast.

Technical Specification :

Digital Display Electronic Universal Testing Machine



Laboratory Instrument India

