



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Digital Sound Sensor

**Product Code :**  
LBNY-0004-780065



**Description :**

Digital Sound Sensor

**Technical Specification :**

Digital Sound Sensor :

Shut-down the automatic off function, sound meter can continuously monitor ambient noise in real time until battery run out of

Mode: quick to know the the highest and lowest sound level around you.

Automatic update the max and min level when monitored the higher or lower sound value.

Digital display (5 digits)

Frequency response from 20 to 7200 Hz, at least

Two stage amplification conditions low level signals

Sound level range from 45 dB to 100 dB

Includes extension cable



**Laboratory Instrument India**

