



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Computer Controlled Electronic Universal Testing Machine

**Product Code :**  
LBNY-0005-1470004



**Description :**

Capability of mechanics testing on various metal materials. e.g. Metal tension, compression, bend, shear, etc.  
Capability of mechanics testing on various nonmetal materials. e.g. Wire tension, plastic bending, cable peeling, etc.  
Closed loop control of channels such as load, displacement, deformation, etc. Smooth switching in between.  
Domestically pioneered process-oriented design of testing control. Provides users with intuitive testing assignment.  
Support for multi-language user-interface extensibility, pre-installed.  
Domestically pioneered open-testing-method settings. Engineers may set up testing projects then operators may employ the projects for testing, greatly reduced the difficulty of the operations and the requirement in the skills and trainings of the operator.  
Dynamic display of load, displacement, deformation and real-time testing curve on the screen.  
Support for displaying a single curve of a specimen or the overlay of multiple curves.

**Technical Specification :**

Computer Controlled Electronic Universal Testing Machine



**Laboratory Instrument India**

