



Email : sales@laboratoryinstrumentindia.com

Product Name :
Microcomputer Control Compression Testing Machine Dual Test

Product Code :
LBNY-0005-1480003



Description :

Domestically pioneered open-testing-method settings. Engineers may set up testing projects then operators may employ the projects for testing, greatly reduced the difficulty of the operations and the requirement in the skills and trainings of the operator.

Stainless steel table provides more convenience in usage and maintenance.

Servo hydraulic controlling system features precision in control, little noise and reduced low energy consumption. Dynamic display of load, displacement, time and real-time testing curve on the screen.

Constant loading speed test supported.

Servo hydraulic system consumes no energy under non-testing states, which saves energy and emits low heat.

Reduced footprint requires little lab space.

Capable of holding on maximum displacement.

Domestically pioneered process-oriented design of testing control. Provides users with intuitive testing control procedures and a wide range of instruction types. Also supports dynamic instruction parameter assignment.

Technical Specification :

Microcomputer Control Compression Testing Machine Dual Test

Laboratory Instrument India

