



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Microcomputer Control Compression Testing Machine

**Product Code :**  
LBNY-0005-1480004



**Description :**

Dynamic display of load, displacement, time and real-time testing curve on the screen.  
Domestically pioneered open-testing-method settings. Engineers may set up testing projects then operators may employ the projects for testing, greatly reduced the difficulty of the operations and the requirement in the skills and trainings of the operator.  
Domestically pioneered process-oriented design of testing control. Provides users with intuitive testing control procedures and a wide range of instruction types. Also supports dynamic instruction parameter assignment.  
Capable of holding on maximum load or displacement.  
Support for compression tests on various metal or non-metal materials to inspect their mechanics performances.  
Closed loop control of channels such as load, displacement, etc. Smooth switching in between.  
Constant loading speed test supported.

**Technical Specification :**

Microcomputer Control Compression Testing Machine



**Laboratory Instrument India**

