



Email : sales@laboratoryinstrumentindia.com

Product Name :
Automatic Compression Testing Machine

Product Code :
LBNY-0005-1480005



Description :

Domestically pioneered process-oriented design of testing control. Provides users with intuitive testing control procedures and a wide range of instruction types. Also supports dynamic instruction parameter assignment. Constant loading speed test supported. Capable of holding on maximum load or displacement. Support for compression tests on various metal or non-metal materials to inspect their mechanics performances. Dynamic display of load, displacement, time and real-time testing curve on the screen. Domestically pioneered highly efficient three-level operation model of testing data processing. Supports open testing data processing, appending custom testing data items and re-inspection of processed testing data. Closed loop control of channels such as load, displacement, etc. Smooth switching in between. Domestically pioneered open-testing-method settings. Engineers may set up testing projects then operators may employ the projects for testing, greatly reduced the difficulty of the operations and the requirement in the skills and trainings of the operator.

Technical Specification :

Automatic Compression Testing Machine



Laboratory Instrument India

