



Email : sales@laboratoryinstrumentindia.com

Product Name :
Microcomputer Control Automatic Compression Testing Machine

Product Code :
LBNY-0005-1480006



Description :

Domestically pioneered open-testing-method settings. Engineers may set up testing projects then operators may employ the projects for testing, greatly reduced the difficulty of the operations and the requirement in the skills and trainings of the operator.

Dynamic display of load, displacement, time and real-time testing curve on the screen.

Constant loading speed test supported.

Capable of holding on maximum load or displacement.

Closed loop control of channels such as load, displacement, etc. Smooth switching in between.

Support for compression tests on various metal or non-metal materials to inspect their mechanics performances.

Technical Specification :

Microcomputer Control Automatic Compression Testing Machine



Laboratory Instrument India