



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Computerized Axial Torque Load Joint Testing Machines

**Product Code :**  
LBNY-0005-1520003



**Description :**

**Features:**

- Automatically sustains when it reaches the specified value.
- Capable of sustaining at the maximum torque, rotation and grip load.
- Dynamic display of torque, rotation and real-time testing curve on the screen.
- Capable of re-analyzing test data and curves, and partial zooming of testing curves at will.
- Integrated with travel limit protection and full-scale overload protection.
- Control over loading speed of torque and rotation, no stepping for loading speed.
- Support for data searching. Returns the corresponding data at y-axis when entered the coordinate at x-axis.
- Support Excel report templates and customized editing of test reports.
- Support for setting a target value on torque, rotation and grip load at will.
- Options for displaying curves including torque-rotation, torque-grip load, grip load-rotation.
- Support saving, instant retrieving and inspection of testing results.

**Technical Specification :**

Computerized Axial Torque Load Joint Testing Machines



**Laboratory Instrument India**

