



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Computerized Axial Torque Load Joint Torsion Testing  
Machine With High Precision Unit

**Product Code :**  
LBNY-0005-1520004



**Description :**

Automatically sustains when it reaches the specified value.  
Capable of sustaining at the maximum torque, rotation and grip load.  
Options for displaying curves including torque-rotation, torque-grip load, grip load-rotation.  
Dynamic display of torque, rotation and real-time testing curve on the screen.  
Control over loading speed of torque and rotation, no stepping for loading speed.  
Integrated with travel limit protection and full-scale overload protection.  
Support saving, instant retrieving and inspection of testing results.  
Capable of re-analyzing test data and curves, and partial zooming of testing curves at will.  
Support for setting a target value on torque, rotation and grip load at will.  
Support for data searching. Returns the corresponding data at y-axis when entered the coordinate at x-axis.  
Support Excel report templates and customized editing of test reports.

**Technical Specification :**

Computerized Axial Torque Load Joint Torsion Testing Machine With High Precision Unit



**Laboratory Instrument India**

