



Email : sales@laboratoryinstrumentindia.com

Product Name :
Computerized Axial Torque Load Joint Torsion Testing
Machine With High Precision Unit

Product Code :
LBNY-0005-1520004



Description :

Automatically sustains when it reaches the specified value.
Capable of sustaining at the maximum torque, rotation and grip load.
Options for displaying curves including torque-rotation, torque-grip load, grip load-rotation.
Dynamic display of torque, rotation and real-time testing curve on the screen.
Control over loading speed of torque and rotation, no stepping for loading speed.
Integrated with travel limit protection and full-scale overload protection.
Support saving, instant retrieving and inspection of testing results.
Capable of re-analyzing test data and curves, and partial zooming of testing curves at will.
Support for setting a target value on torque, rotation and grip load at will.
Support for data searching. Returns the corresponding data at y-axis when entered the coordinate at x-axis.
Support Excel report templates and customized editing of test reports.

Technical Specification :

Computerized Axial Torque Load Joint Torsion Testing Machine With High Precision Unit



Laboratory Instrument India

