



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Computerized Electronic Wire Torsion Testing Machine

**Product Code :**  
LBNY-0005-1520005



**Description :**

**Features:**

Automatic stop when reaches target rotation turns or specimen breaks.  
Distance between grip jaw continuously adjustable by moving the ball screw.  
Stably and reliably controlled through human interface and PLC.  
Stainless table makes it easier to clean.  
TFT true-color LED screen to display current rotation, target rotation turn and configured speed.  
No stepping in speed adjustment. Automatic data recording. Rotation speed and turns configurable.  
Equipped with translucent safety shield to prevent specimen from popping out when it breaks.  
More reliable hydraulic gripping means for wires in large diameter.

**Technical Specification :**

Computerized Electronic Wire Torsion Testing Machine



**Laboratory Instrument India**