



Email : sales@laboratoryinstrumentindia.com

Product Name :
Computer Control High Precision Torsion Testing Machine

Product Code :
LBNY-0005-1520007



Description :

Dynamic display of torque, rotation and real-time testing curve on the screen.
Automatically sustains when it reaches the specified value.
Capable of sustaining at the maximum torque, rotation.
Control over loading speed of torque and rotation, no stepping for loading speed.
Support saving, instant retrieving and inspection of testing results.
Integrated with travel limit protection and full-scale overload protection.
Capable of re-analyzing test data and curves, and partial zooming of testing curves at will.
Support for setting a target value on torque, rotation and grip load at will.
Options for displaying curves including torque-rotation.
Support for data searching. Returns the corresponding data at y-axis when entered the coordinate at x-axis.
Support Excel report templates and customized editing of test reports.

Technical Specification :

Computer Control High Precision Torsion Testing Machine



Laboratory Instrument India

