



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital Display Torsion Tester

Product Code :
LBNY-0005-1520008



Description :

TFT true-color LED screen to display current rotation, target rotation turn and configured speed.
No stepping in speed adjustment. Automatic data recording. Rotation speed and turns configurable.
Equipped with translucent safety shield to prevent specimen from popping out when it breaks.
More reliable hydraulic gripping means for wires in large diameter.
Automatic stop when reaches target rotation turns or specimen breaks.
Distance between grip jaw continuously adjustable by moving the ball screw.
Stably and reliably controlled through human interface and PLC.
Stainless table makes it easier to clean.

Technical Specification :

Digital Display Torsion Tester



Laboratory Instrument India