



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Horizontal Bend Testing Equipment

**Product Code :**  
LBNY-0005-1530003



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Main applications: quality control department, iron forging, civil aerospace, higher education, science research department, arbitration department, technical supervision department and etc.

This product features an appropriate structure and layout, as well as an ergonomic design, making it more convenient to operate. This product emits little noise, consumes little energy and fulfills the requirements of energy saving and environment protection. This product adopts the plug-in flexural center installation style, presenting an efficient design that requires a small footprint.

Flexure tests can be performed automatically or manually. The LWC series of programmable flexure testing machine is mainly designed for flexure tests on steel plate specimen. This product features an integrated main frame. The vertical and the horizontal two-way hydraulic cylinders are designed to form a complete unit. The span of the horizontal rollers may be adjusted freely and simultaneously. The two rollers are capable of 180°, ° flexure tests.

**Technical Specification :**

Horizontal Bend Testing Equipment



**Laboratory Instrument India**

