



Email : sales@laboratoryinstrumentindia.com

Product Name :
Programmable Flexure Testing Machine

Product Code :
LBNY-0005-1530004



www.laboratoryinstrumentindia.com

Description :

The vertical and the horizontal two-way hydraulic cylinders are designed to form a complete unit. The span of the horizontal rollers may be adjusted freely and simultaneously. The series of programmable flexure testing machine is mainly designed for flexure tests on steel plate specimen. The two rollers are capable of 180° flexure tests. Flexure tests can be performed automatically or manually. This product features an integrated main frame.

This product emits little noise, consumes little energy and fulfills the requirements of energy saving and environment protection. This product features an appropriate structure and layout, as well as an ergonomic design, making it more convenient to operate. This product adopts the plug-in flexural center installation style, presenting an efficient design that requires a small footprint.

Main applications: Quality control department, iron forging, civil aerospace, higher education, science research department, arbitration department, technical supervision department and etc.

Technical Specification :

Programmable Flexure Testing Machine



Laboratory Instrument India

