



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Computerized Hydraulic Bend Testing Machine

**Product Code :**  
LBNY-0005-1530006



**Description :**

Safe: specimen bending directly into sealed test table, safe to operate.  
Dust-proof: unique dust resistant unit protects ball screw from oxide coating.  
Small footprint saves space for testing. Operating: digital buttons, touch-screen interactions, performing operations through only some simple clicks.  
Operating: digital buttons, touch-screen interactions, performing operations through only some simple clicks.  
Effortless: requiring only one operator, who clicks the buttons to perform tests.  
Effortless: requiring only one operator, who clicks the buttons to perform tests.

**Technical Specification :**

Computerized Hydraulic Bend Testing Machine



**Laboratory Instrument India**