



Email : sales@laboratoryinstrumentindia.com

Product Name :
Computerized Hydraulic Bend Testing Machine

Product Code :
LBNY-0005-1530006



Description :

Safe: specimen bending directly into sealed test table, safe to operate.
Dust-proof: unique dust resistant unit protects ball screw from oxide coating.
Small footprint saves space for testing. Operating: digital buttons, touch-screen interactions, performing operations through only some simple clicks.
Operating: digital buttons, touch-screen interactions, performing operations through only some simple clicks.
Effortless: requiring only one operator, who clicks the buttons to perform tests.
Effortless: requiring only one operator, who clicks the buttons to perform tests.

Technical Specification :

Computerized Hydraulic Bend Testing Machine



Laboratory Instrument India