



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electro Hydraulic Servo Microcomputer Control Static Load
Anchor Testing Machine

Product Code :
LBNY-0005-1550006



www.laboratoryinstrumentindia.com

Description :

Features:

Dynamic display of load, displacement, time and real-time testing curve on the screen.
Closed loop control of channels such as load, displacement, etc. Smooth switching in between.
Servo-controlled loading. Controllable pre stress. Load is applied to specimen uniformly.
Support for holding at constant load or constant displacement.
Software online notification of current working status and operating procedures.
Capable of loading at constant speeds of various levels.
Capable of automatic loading or holding by loading program settings.
Domestically pioneered process-oriented design of testing control. Provides users with intuitive testing control procedures and a wide range of instruction types. Also supports dynamic instruction parameter assignment.
Engineers may set up testing projects then operators may employ the projects for testing, greatly reduced the difficulty of the operations and the requirement in the skills and trainings of the operator.

Technical Specification :

Electro Hydraulic Servo Microcomputer Control Static Load Anchor Testing Machine

Laboratory Instrument India

