



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Computerized Hydraulic Servo Static Load Anchor Testing Machine

**Product Code :**  
LBNY-0005-1550007



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

**Features:**

Servo-controlled loading. Controllable pre stress. Load is applied to specimen uniformly.  
Dynamic display of load, displacement, time and real-time testing curve on the screen.  
Closed loop control of channels such as load, displacement, etc. Smooth switching in between.  
Domestically pioneered process-oriented design of testing control. Provides users with intuitive testing control procedures and a wide range of instruction types. Also supports dynamic instruction parameter assignment.  
Capable of loading at constant speeds of various levels.  
Capable of automatic loading or holding by loading program settings.  
Engineers may set up testing projects then operators may employ the projects for testing, greatly reduced the difficulty of the operations and the requirement in the skills and trainings of the operator.  
Support for holding at constant load or constant displacement.  
Software online notification of current working status and operating procedures.

**Technical Specification :**

Computerized Hydraulic Servo Static Load Anchor Testing Machine

---

# Laboratory Instrument India

