



Email : sales@laboratoryinstrumentindia.com

Product Name :
Optical Cable Wire Testing machine with Micrcomputer Control

Product Code :
LBNY-0005-1570007



www.laboratoryinstrumentindia.com

Description :

This product features framed feed-through structure, making it easier to install optical cable specimen. Integrated non-contact transducer detects torsion rotation, avoiding the abrasion caused by contact transducers, ensuring precision controls mean while providing a long operating life. Remote communication ports are installed to make it easier for remote controlling operation on a computer.

Technical Specification :

Torsion length: 1000mm.
Torsion times: 19999 times, configured at will.
Torsion speed: 530 times/min, configured at will.
Torsion rotation: 90,180,360 for choice.
Total hammer weight: 22.5kg. Weights: 5kg4,2.5kg1
Hammer bracket weight: 5kg.



Laboratory Instrument India

