



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wireless Force and Acceleration Sensor

Product Code :
LBNY-0004-780069



Description :

Wireless Force and Acceleration Sensor

Technical Specification :

Wireless Force and acceleration Sensor
Force range: -50 to +50 N with a 0.03N resolution
Tare button to zero out small forces
3-axis accelerometer
3-axis gyroscope
Acceleration range: -16 to 16 g.
Supplied with adapter to connect sensor to PC or laptop.



Laboratory Instrument India