



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pitot Apparatus

Product Code :
LBNY-0005-1580005



Description :

The Pitot Apparatus consists of a total head tube and a static head tube for the Pitot, clear acrylic pipe, inlet and outlet, and an inclined manometer. Pitot Apparatus is for a study of velocity distribution across pipe section for a water flow in pipe, and for measurement of flow velocity. The total head tube can be moved up or down vertically with a scale for indication of the tube position. Differential pressure is measured by a multiple inclination manometer. Pitot Apparatus is to be used with Hydraulics Bench.

Technical Specification :

Acrylic tube: 40 mm. diameter.
Inclined water manometer: 0-450.
Total head probe: 1 ea.



Laboratory Instrument India