



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydrostatics Panel

Product Code :
LBNY-0005-1580006



Description :

The unit can be placed on a suitable bench top and requires outside water supply to be stored in the provided storage tank. The Hydrostatics Panel equipment is a bench top unit for the study of basic properties of fluids in hydrostatic conditions. Components are placed on the panel.

Technical Specification :

The Hydrostatics Panel equipment on the panel comprises:

Barometer: Hygrometer.

Thermometer: Hydrometer

Inclined water manometer.

Specific gravity bottle: Triple beam balance

Beaker, graduated cylinder and pipette: Archimedes apparatus.

Stainless steel vernier hook gauge: Pascals apparatus.

Capillary viscometer: Falling sphere viscometer

Archimedes Principle Apparatus: Storage tank

Capillary attraction plates : Capillary tubes.



Laboratory Instrument India

