



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hydrostatics Set

**Product Code :**  
LBNY-0005-15800011



#### **Description :**

Hydrostatics Set is a reduced version of the hydrostatics bench designed to be used with the Hydraulics Bench.

Experimental Capability of Hydrostatics Set:

Manometry and Hares Tube

Capillarity

Demonstration of Pascals Law.

The unit consists of a backboard unit, mounted directly onto the with:

Hare tube, 600 mm long

Manometer manifold block

2x Fluid manometer, 600 mm long

Pascals apparatus, comprising glass tubes of differing shapes and varying cross-sectional area sealed by a pressure disc actuated by a weighing beam

Capillarity apparatus comprising set of capillary tubes in frame and trough

Capillarity apparatus comprising two glass plates in frame and trough

#### **Technical Specification :**

Dimensions And Weights:

Length: 36cm

Width: 10cm

Height: 75cm

Nett weight: 9kg.



**Laboratory Instrument India**