



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Force Sensor Balance Stand

**Product Code :**  
LBNY-0004-780071



**Description :**

Force Sensor Balance Stand

**Technical Specification :**

Force Sensor Balance stand

Stand used for Pan Balance, measure Buoyant force, demonstrate the change in normal force as the angle of support changes.

Force Sensor to this convenient stand and students have an electronic balance for a wide variety of physics experiments but with the added feature of using the data collection and analysis power of software.

Connect an Acceleration Sensor for studies of angle vs. normal force. The mounting screws and balance pan can be stored on the stand when not in use.



**Laboratory Instrument India**