



Email : sales@laboratoryinstrumentindia.com

Product Name :
Distribution Trainer

Product Code :
LBNY-0005-15900012



Description :

Nett dimensions and weight:

1850 mm long x 1870 mm high x 960 mm front to back and 670kg

Packed volume and weight: 5 m³ and 800 kg

Three-phase loads:

220 V star-connected resistive, inductive and capacitive

110 V delta-connected resistive

Protection relays:

1 x single phase over current

1 x three phase over current

1 x three-phase over current with voltage

Transformers:

2 x low-voltage single-phase

1 x line power 5 kVA three-phase delta-star reduces the incoming supply down to 220/230 V

1 x line receive 3 kVA three-phase delta-star with starpoint link to earth 220 V/110 V three-phase. Includes tap selection on its primary windings of $\hat{A}\pm 10\%$, $\hat{A}\pm 7.5\%$, $\hat{A}\pm 5\%$, $\hat{A}\pm 2.5\%$, and 0%.

Low voltage supplies:

30 volts d.c. or 0-30 volts a.c.

Technical Specification :

Distribution Trainer



Laboratory Instrument India