



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wireless Conductivity Sensor

Product Code :
LBNY-0004-780072



Description :

Wireless Conductivity Sensor

Technical Specification :

Wireless Conductivity sensor:
Range: 0 to 20000 μ S/cm
Resolution: at least 0.1 μ S/cm
Accuracy: \pm 10% of value from 200 μ S/cm to 20,000 μ S/cm
Below 200 μ S/cm: qualitative
Temperature compensation: 0 to around 80°C
Temperature accuracy: \pm 0.5 C
Includes adapter to connect sensor to laptop or PC



Laboratory Instrument India