



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Wireless Conductivity Sensor

**Product Code :**  
LBNY-0004-780072



**Description :**

Wireless Conductivity Sensor

**Technical Specification :**

Wireless Conductivity sensor:  
Range: 0 to 20000  $\mu$ S/cm  
Resolution: at least 0.1  $\mu$ S/cm  
Accuracy:  $\pm$  10% of value from 200  $\mu$ S/cm to 20,000  $\mu$ S/cm  
Below 200  $\mu$ S/cm: qualitative  
Temperature compensation: 0 to around 80°C  
Temperature accuracy:  $\pm$  0.5 C  
Includes adapter to connect sensor to laptop or PC



**Laboratory Instrument India**