



Email : sales@laboratoryinstrumentindia.com

Product Name :
Variable Three Phase Power Supply

Product Code :
LBNY-0005-15900017



Description :

Power supply unit for variable three phase voltage suitable for supplying AC machines.
Multifunctional digital instrument that supplies voltage, current and power
Ac output: 3 x 0 ... 380 V, 8 A
Start and stop push-buttons
Motor protection circuit breaker: 6.3 to 10 A
16 A, 30 mA differential magneto-thermal main switch
Key operated emergency push-button
Modbus interface connectors.
Supply voltage: three-phase from mains.
The output voltage is set by a rotary knob with 0-100% scale.

Technical Specification :

Variable Three Phase Power Supply



Laboratory Instrument India

