



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Electrical Machines & Drives Program

**Product Code :**  
LBNY-0005-15900019



**Description :**

Accurate torque meter systems. The new Torque Sensor uses the latest technology, without need for regular calibration.  
A handy physical size.  
Possibility to use standard instruments.  
A comprehensive range of test machines with accessories, to show typical characteristics, valid also for larger machines.

**Technical Specification :**

Electrical Machines & Drives Program



**Laboratory Instrument India**