



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Wireless Heart Rate Sensor

**Product Code :**  
LBNY-0004-780073



**Description :**

Wireless Heart Rate Sensor

**Technical Specification :**

Wireless Heart rate Sensor:

With these wireless hand grips, conducting physiology labs on the cardiovascular system or homeostasis is easier than ever before.

Continuously monitor heart rate during exercise, or use the sensor to take initial and final measurements with fast and reliable heart rate detection.

Used to demonstrate the cardiovascular system

Come with hand-grip accessory

Include adapter to connect the sensor to laptop or PC.



**Laboratory Instrument India**