



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Set Of Sectioned Components For Liquid Pumping

**Product Code :**

LBNY-0005-1600007



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Special explanatory technical plates show each sectioned component explaining its operation with notes, drawings and technical data. The Set Of Sectioned Components For Liquid Pumping includes some of the most important components, used in water works, for fluid pumping, that is pumps of different types. The use of contrasting colours according to a pre-established code simplifies the informative process as well as the learning one. In sectionizing each component the best care has been taken to evidence all the interesting constructive detailed features.

**Technical Specification :**

Set Of Sectioned Components For Liquid Pumping



**Laboratory Instrument India**